

# CORONAVIRUS DISEASE 2019 ANTIBODY (IGG/IGM)

COMBINED TEST KIT DIXION



**DIXION**

# CORONAVIRUS DISEASE 2019 ANTIBODY (IGG/IGM)

- ❖ RAPID Results within 5-10 minutes
- ❖ EASY Only 4 steps taken


- ❖ DIVERSE Serum, Plasma, Whole blood
- ❖ RELIABLE IgG/IgM tested at one time

## Specification:

Product name	Coronavirus Disease 2019 Antibody (IgM/IgG) Combined Test Kit DIXION
Test principle	Lateral flow immunoassay
Sample	Whole blood, Serum, Plasma
Time to report	5-10 minutes
Operation temperature	Room temperature
Pack	25 tests
Material provided	25 individual sealed pouches, 25 disposable droppers

DIXION IgM/IgG antibody combined rapid test for COVID-19 is used to qualitatively detect IgG and IgM antibodies of the novel coronavirus (SARS-CoV-2) in human serum, plasma or whole blood in vitro. These biomarkers indicate that a person's immune system has responded to infection with the COVID-19 virus.

This innovative rapid test kit is developed with lateral flow immunoassays technology in high sensitivities and specificities. It is quite easy to use and ideal for screening check in epidemics control of COVID-19 virus.



Allow the test cassette, specimen, buffer to equilibrate to room temperature (15-30°C) prior to testing.

- Remove the test cassette from the sealed foil pouch and use it as soon as possible. Best results will be obtained if the assay is performed immediately after opening the foil pouch.
- Collect blood sample, transfer 10 µL of serum / plasma specimen or 20 µL (1 drop) of whole blood specimen to the sample well.
- Add about 80µL (about 2-3 drops) of buffer and start the timer.
- Wait for the colored line(s) to appear. Read the result at 5-10 minutes.
- Do not interpret the result after 15 minutes.

# DIXION

Vertrieb medizinischer Geräte GmbH

Couvenstr. 6, 40211 Düsseldorf, Germany

Tel. +49 (0)211 38838868

Fax +49 (0)211 38838697

[www.dixion.de](http://www.dixion.de)  
[info@dixion.de](mailto:info@dixion.de)